

Phase Resolved Doppler Optical Coherence Tomography (prDOCT) Data Processing Guide

Jianbo Tang

Aug 7, 2019

Introduction:

This guide is for post data processing of prDOCT which outputs 3D vascular blood flow velocity.

Please cite the following references^{1,2}:

1. Tang, J., Erdener, S. E., Fu, B. & Boas, D. A. Capillary red blood cell velocimetry by phase-resolved optical coherence tomography. *Opt. Lett.* **42**, 3976–3979 (2017).
2. Tang, J., Erdener, S. E., Sunil, S. & Boas, D. A. Normalized field autocorrelation function-based optical coherence tomography three-dimensional angiography. *J. Biomed. Opt.* **24**, 036005 (2019).

I. Data acquisition

OCT-based M-mode data acquisition, i.e. repeat Ascan at each X-Y scanning location for a certain period. The data should be saved sequentially as a 1D array (ASCII int16) and named as: RAW-nk-nxRpt-nx-nyRpt-ny-iC, e.g. RAW-1024-100-00400-001-400-1.dat.

II. Input

```
% 1D array spectrum, nK*nXrpt*nX*nY, data format: ASCII int16
% nK: spectrum pixel (camera elements); nXrpt: Ascan repeat;
% nX: number of Ascans per Bscan; nY: number of Bscans for the whole volum
% NOTE: the raw data for the whole volume is usually very large, it's
recommended to process chunk by chunk
% PRSinfo: processing information
% PRSinfo.FWHM: Full width at Half Maxim, Amplitude, [transverse, axial], m
% PRSinfo.fAline: DAQ Aline rate, Hz
% PRSinfo.Lam: [light source center, wavelength bandwidth], m
% PRSinfo.Dim: [nz,nx,nyPchk,nTau]
% PRSinfo.intDk: OCT lambda to k space interpolation factor (calibration is
required)
```

Example data:

<https://drive.google.com/open?id=168HD4IKt0K97g09zus6H9h7lAyO0jOBZ>
https://drive.google.com/open?id=1QvTO_41cPN3_wM9wxCh9NECv_hypVZPC

III. Output

```
% Vz, [nz,nx,ny]
```

I. CPU calculation-based sub-functions

```
% subFunctions:
% function [Dim, fNameBase, fIndex]=GetNameInfoRaw(filename0)
% function DAT= ReadDat_int16(filePath, Dim, iseg, ARpt_extract,RptBscan)
% function RR = DAT2RR(Dat, intpDk)
% function GG = RR2g1(RR, PRSinfo)
% function [Vz, aRBC]=RR2Vz(RR, PRSinfo)
```


II. GPU calculation-based sub-functions

Note: the minimal GPU memory requirement is 16 GB.

```
% subFunctions:
% function [Dim, fNameBase, fIndex]=GetNameInfoRaw(filename0)
% function DAT= ReadDat_int16(filePath, Dim, iseg, ARpt_extract,RptBscan)
% function RR = DAT2RR_GPU(Dat, intpDk)
% function GG = RR2g1_GPU(RR, PRSinfo)
% function [Vz, aRBC]=RR2Vz_GPU(RR, PRSinfo)
```


III. Main_g1OCTA data processing

III.1 select file



III.2 data processing parameter

Specify the number of chunks for data processing (split large data size into small chunks), and the g1 calculation parameters (nTau and nt). Double check the intDk and Aline Rate.




A MATLAB dialog box with the following fields and values:


- Num Ysegs (nY=400, nX=400, NxRpt=150, NyRpt=1): 40
- RptA_Start (nARpt Process): 1
- RptA_Interval (nARpt Process): 1
- RptA_n (nARpt Process): 150
- N_Aline(DAQ): 400
- N_Bscan(DAQ): 400
- N_xRpt(DAQ): 150
- N_yRpt(DAQ): 1
- intDk: -0.43
- Aline rate (kHz): 46

Buttons: OK, Cancel

III.3 select the axial data processing range [surface, zStart, zEnd]



III.4 g1OCTA result plot (GUI_PLOT_DLSOCT)



1. Load the saved prVz data
2. Use the slider or zStart+zStack to check single or MIP (maxim intensity projection) en face plane. Cyan color means descending flow
3. Use 'Plot MIP&SideView' to plot a MIP for certain depth range (set SideView(N:0; Y:1) to 0). Or plot XY, YZ, XZ, and MIP figures by set SideView(N:0; Y:1) to 1.

